

Automatic Content Progression

Previously, it was possible to schedule Case Collections within a session to automatically become available to users based on time.

New in PathoTrainer 2.2.0 is the ability to set episodic content release also by progression, ensuring that your participants don't rush ahead or skip out crucial training & assessment exercises.

The video below will guide you through how to set this up for yourself

[automatic_progression.mp4](#)

From:

<https://docs.pathomation.com/pathotrainer/> - **Pathotrainer**

Permanent link:

https://docs.pathomation.com/pathotrainer/doku.php?id=automatic_content_progression

Last update: **2025/11/25 18:24**

