


# Adding Cases and Case Collections to a Project

Case collections are how you organize your content to be shown in courses.

## Creating a case collection in a project


You can create a Case Collection via the [slides and images](#) option in the sidebar of Pathotrainer and assign that case collection to a project during setup.

Or you can directly add a case collection to a project from within the project details window (shown below). If you have not added any case collections to your project, you will see the following bluetext (1) prompting you to do so:




The screenshot shows the Pathotrainer interface for a project named 'Pathoproject\_T3'. The sidebar on the left contains navigation options: Home, Users, Projects, Slides & Images, Quizzes, Training sessions, Forms, System, Backend, and Documentation. The main content area displays the project details for 'Pathoproject\_T3'. The details are organized into sections: Title (Pathoproject\_T3), Logo, State (Active), Description, Contact info, Viewer theme (Default), Project Owners, Case collections, and Training sessions. The Case collections section shows the message 'There are no case collections in this project' followed by a blue link 'Setup a case collection'. A red arrow points to this link, and a red '1' is placed next to it, indicating the prompt mentioned in the text.


If you have already added some case collections, you can click on the bluetext name of the case collection(s) (1) to view the case collection:



then click Case Collections (2),




and then Create (3).



### Creating a New Case Collection

After hitting the “Create” button, you will be presented with the option to define some details about your new collection:


1. Uploading a logo image
2. Assigning the case to a project
3. Giving the case a name
4. Assigning a scoring form for your trainees
5. Assigning a clinical information form for your trainees
6. Options to toggle slide label and measurement visibility
7. Additional comment templates



### Adding Cases


Now that you've created a case, its time to add some cases.

To do this, navigate to the “Cases” tab and click the “Create”



You will find yourself in the Editing Cases tab. Here you can set up various case properties:

- A name for the case
- What slide information will be displayed: no labels, filenames, the caption
- Comments
- Select your slide(s) from your PMA core inventory
- Save your slides



### A Note on Adding Content

You can use any slides from your [slide inventory](#) in Pathotrainer case collections. The slide inventory is managed by PMA.core - for guides on how to import content there, you might find these links useful:

- [https://docs.pathomation.com/pma.core/3.0.2/doku.php?id=minimal\\_user\\_management](https://docs.pathomation.com/pma.core/3.0.2/doku.php?id=minimal_user_management)
- <https://docs.pathomation.com/pma.core/3.0.2/doku.php?id=upload>

From:  
<https://docs.pathomation.com/pathotrainer/> - **Pathotrainer**

Permanent link:  
[https://docs.pathomation.com/pathotrainer/doku.php?id=cases\\_and\\_case\\_collections&rev=1701688161](https://docs.pathomation.com/pathotrainer/doku.php?id=cases_and_case_collections&rev=1701688161)

Last update: **2023/12/04 14:09**

