



# Edit A Project

After hitting “Create” in the [Create a Project](#) menu, or navigating to an existing project by clicking on a project's name (1) or clicking “Edit” from the dropdown menu, you will be taken to the editing page.



Here you have access to several tabs which allow you to:

1. Ordered List Itemsee an overview of the project Details
2. Edit the project options chosen during [Create a project](#)
3. Add additional material
4. Add links
5. Perform integrity checks
6. Create custom annotations



Project Details

Here you can get an overview of what content and training sessions is already in the project. If no content has been added yet, you are made aware by the orange text and can click the bluetext link to get started

Patho\_Project

[Home](#) / [Projects](#) / [Patho\\_Project](#)

Details

Edit

Additional material

Links

Integrity check

Annotations toolbar

Title

Patho\_Project

Logo

State

Active

Description

A project to show how easy it is to start a project!

Contact info

Viewer theme

Modern

Project Owners

Case collections

There are no case collections in this project

[Setup a case collection](#)

Training sessions

There are no training sessions in this project

[Setup a training session](#)

Allowed interaction modes ?

☒ Board

☒ Browse

☒ Clinical information edit

☒ Consensus score edit

☒ Consensus view

☒ Hidden

☒ Locked

☒ Review

☒ Self review

☒ Self test

☒ Test active

Back to list

Archive

From:

<https://docs.pathomation.com/pathotrainer/> - **Pathotrainer**

Permanent link:

[https://docs.pathomation.com/pathotrainer/doku.php?id=edit\\_a\\_project&rev=1701348695](https://docs.pathomation.com/pathotrainer/doku.php?id=edit_a_project&rev=1701348695)

Last update: **2023/11/30 15:51**

