



# Edit A Project

After hitting “Create” in the [Create a Project](#) menu, or navigating to an existing project by clicking on a project's name (1) or clicking “Edit” from the dropdown menu, you will be taken to the editing page.



Here you have access to several tabs which allow you to:

1. Ordered List Itemsee an overview of the project Details
2. Edit the project options chosen during [Create a project](#)
3. Add additional material
4. Add links
5. Perform [integrity checks](#)
6. Create [custom annotation tools](#)



## Project Details

Here you can get an overview of what content and training sessions is already in the project. If no content has been added yet, you are made aware by the orange text and can click the bluetext link to get started

Patho\_Project

[Home](#) / [Projects](#) / [Patho\\_Project](#)

Details	Edit	Additional material	Links	Integrity check	Annotations toolbar
---------	------	---------------------	-------	-----------------	---------------------

Title

Patho\_Project

Logo

State

Active

Description

A project to show how easy it is to start a project!

Contact info

Viewer theme

Modern

Project Owners

Case collections

There are no case collections in this project

[Setup a case collection](#)

Training sessions

There are no training sessions in this project

[Setup a training session](#)

Allowed interaction modes ?

☒ Board

☒ Browse

☒ Clinical information edit

☒ Consensus score edit

☒ Consensus view

☒ Hidden

☒ Locked

☒ Review

☒ Self review

☒ Self test

☒ Test active

Back to list

Archive

## Adding Material and Links

In these tabs you can add material and links, and see anything that has already been added

Additional material for Patho\_Project

Home / Projects / Patho\_Project / Additional material

Details

Edit

Additional material

Links

Integrity check

+ Add

Recycle bin

There are no attachments

Links of Patho\_Project

Home / Projects / Patho\_Project / Links

Details

Edit

Additional material

Links

Integrity check

+ Add

Filter Patho\_Project

There are no links

From:  
<https://docs.pathomation.com/pathotrainer/> - **Pathotrainer**

Permanent link:  
[https://docs.pathomation.com/pathotrainer/doku.php?id=edit\\_a\\_project&rev=1701353231](https://docs.pathomation.com/pathotrainer/doku.php?id=edit_a_project&rev=1701353231)

Last update: **2023/11/30 17:07**

