

Home page

Upon successful log in you arrive on the Home Page.

Here you can see: In this page you can select any of the available sessions to join. If you are registered to only one session, the system will automatically redirect you to the particular session's page.

1. A list of Projects in which you have been registered
2. The names of sessions you have been registered in and your role in them
3. How many cases you have scored/evaluated in this session
4. The date(s) that the course is Active. If the course is currently active, the text will be green. If you are outside the active dates, the text will display red.

You can enter a training session by clicking the 'Join' button

The screenshot shows the 'Home' page of the Pathotrainer application. At the top, there is a header bar with 'PMA.control test environment' on the left and a user profile icon labeled 'User' on the right. Below the header, the page title is 'Home'. The main content area features a welcome message: 'Welcome to the HER2 IHC Interpretation training class'. This is followed by two paragraphs of text describing the session's focus on HER2 immunohistochemistry (IHC) interpretation and its correlation with other diagnostic techniques like FISH. Below the text, there is a section titled 'Patho_Project' which contains a card for a 'Training Session Test (Trainee)'. This card displays '0/16 cases evaluated', the dates 'Wed. 3 Jan 2024 - Thu. 4 Jan 2024', and the lead 'Lead By Sea Bow'. A 'Join >' button is located at the bottom right of the card. Red arrows with numbers 1 through 4 point to specific elements: 1 points to the 'Patho_Project' title, 2 points to the 'Training Session Test (Trainee)' title, 3 points to the '0/16 cases evaluated' text, and 4 points to the date range. Below the card, the text 'HER2 IHC' is visible, followed by a logo for 'HER2 PRECISION MEDICINE'.

From:

<https://docs.pathomation.com/pathotrainer/> - Pathotrainer

Permanent link:

https://docs.pathomation.com/pathotrainer/doku.php?id=home_page

Last update: **2024/01/03 15:48**

