


## Home page

Upon successful log in, a list of sessions in which you have been registered as a participant by the project owners or administrators of the platform, is presented. In this page you can select any of the available sessions to join. If you are registered to only one session, the system will automatically redirect you to the particular session's main page.



From:  
<https://docs.pathomation.com/pathotrainer/> - **Pathotrainer**

Permanent link:  
[https://docs.pathomation.com/pathotrainer/doku.php?id=home\\_page&rev=1596716340](https://docs.pathomation.com/pathotrainer/doku.php?id=home_page&rev=1596716340)

Last update: **2020/08/06 15:19**

