

Homepage (Supervisor)

Upon successful log in, you will see next to the list of sessions in which you have been registered as a user as well as the sessions in which you are registered as **supervisor**.

In this page you can see Projects with sessions you can access (1)

See Training Session names that you are able to be a participant of (2)

See the period when the training session will be live and accessible (3) - this text will be green when the session is active, and red when inactive.

PMA.control test environment v2.0.0.1336

Home

Projects

Slides & Images

Training sessions

Forms

Documentation

Welcome to the HER2 IHC Interpretation training class

In this session, we will dive into the intricate world of HER2 Immunohistochemistry (IHC) interpretation, an epidermal growth factor receptor, is a crucial biomarker that aids in determining the prognosis and treatment testing, including the interpretation of staining patterns, intensity, and scoring systems.

We will also discuss the challenges and pitfalls associated with HER2 IHC interpretation and examine the co situ hybridization (FISH). By the end of this training, you will have a comprehensive understanding of HER2 I practice.

Patho_Project

Training Session Test (Supervisor)

Wed, 3 Jan 2024 - Thu, 4 Jan 2024

Lead By Sea Bow

Join >

From:

<https://docs.pathomation.com/pathotrainer/> - Pathotrainer

Permanent link:

https://docs.pathomation.com/pathotrainer/doku.php?id=home_page_supervisor

Last update: **2024/01/03 13:51**

