

# Homepage

Upon successful log in, you will see next to the list of sessions in which you have been registered as a participant/trainee as well as the sessions in which you are registered as supervisor (trainer).

In this page you can select any of the available sessions to join.

## Home

You have been registered in the following training sessions. Please choose one to join.

Filter

All projects

Any role

DEMO KIDN		
<b>Demo session 1</b> (Trainee) 0/3 cases evaluated	<a href="#">Join</a>	----- Star
<b>Demo session 2</b> (Supervisor)	<a href="#">Join</a>	Lead by wim waelput

From:

<https://docs.pathomation.com/pathotrainer/> - **Pathotrainer**

Permanent link:

[https://docs.pathomation.com/pathotrainer/doku.php?id=home\\_page\\_supervisor&rev=1601048536](https://docs.pathomation.com/pathotrainer/doku.php?id=home_page_supervisor&rev=1601048536)

Last update: **2020/09/25 18:42**

