

Homepage

Upon successful log in, you will see next to the list of sessions in which you have been registered as a user as well as the sessions in which you are registered as **supervisor** (trainer).

In this page you can select any of the available sessions to join.

The screenshot shows the 'Home' page of the Pathotrainer application. At the top, it says 'Home'. Below that, a message reads: 'You have been registered in the following training sessions. Please choose one to join.' There are two filter dropdowns: 'Filter All projects' and 'Any role'. Below the filters is a table of training sessions. The table has a dark header with the text 'DEMO KIDN'. The first row is 'Demo session 1 (Trainee)' with '0/3 cases evaluated' and a 'join' button. The second row is 'Demo session 2 (Supervisor)' with 'Lead by wim waelput' and a 'join' button. A red arrow points to the 'Demo session 2 (Supervisor)' row.

DEMO KIDN	
Demo session 1 (Trainee) 0/3 cases evaluated	join
Demo session 2 (Supervisor)	join

From:
<https://docs.pathomation.com/pathotrainer/> - Pathotrainer

Permanent link:
https://docs.pathomation.com/pathotrainer/doku.php?id=home_page_supervisor&rev=1635363241

Last update: 2021/10/27 22:34

