

# Live mode

From time to time the trainer/supervisor of a session can take over control of your the viewport of all participants of the sessions. The trainer/supervisor will navigate through the cases and individual whole slide images.

Live session in action video:

(under construction)

The Live mode is supported by a separate Pathomation software component called [PMA.live](#).

From:

<https://docs.pathomation.com/pathotrainer/> - **Pathotrainer**

Permanent link:

[https://docs.pathomation.com/pathotrainer/doku.php?id=live\\_mode&rev=1647959693](https://docs.pathomation.com/pathotrainer/doku.php?id=live_mode&rev=1647959693)

Last update: **2022/03/22 17:34**

