

# Quick Start

This quick-start guide is intended to point you in the right direction to create a basic virtual classroom environment in PathoTrainer in which your Trainees will be able to view cases and input data.

This guide assumes that you have an administrator level account.

1. [Upload slides](#)
2. [Create a Project](#)
3. [Create a Case Collection](#)
4. [Create a Training Session](#)
5. [Create Users](#)
6. Optional: [Create](#) and [Assign](#) a form if you want to collect scoring data from your trainees
7. [Invite participants!](#)

From:

<https://docs.pathomation.com/pathotrainer/> - **Pathotrainer**

Permanent link:

[https://docs.pathomation.com/pathotrainer/doku.php?id=quick-start\\_guide](https://docs.pathomation.com/pathotrainer/doku.php?id=quick-start_guide)

Last update: **2025/11/27 11:12**

