

# Quick Start

This quick-start guide is intended to point you in the right direction to create a basic virtual classroom environment in PathoTrainer in which your Trainees will be able to view cases and input data.

This guide assumes that you have an administrator level account.

1. Upload slides
2. Create a Project
3. Create a Case Collection
4. Create a Training Session
5. Create Users
6. Optional: Create a Form if you want to collect scoring data from your trainees
7. Invite your users!

From:

<https://docs.pathomation.com/pathotrainer/> - **Pathotrainer**

Permanent link:

[https://docs.pathomation.com/pathotrainer/doku.php?id=quick-start\\_guide&rev=1763130159](https://docs.pathomation.com/pathotrainer/doku.php?id=quick-start_guide&rev=1763130159)

Last update: **2025/11/14 17:22**

