

# Regions

In the Regions page of the System section of the administrators toolbar, you will see your servers and geo-replication instances.

In this table you can see

1. The name of your instance
2. Where it is located
3. The web address of your PMA.core
4. The version and license information
5. The web address of your PMA.live
6. The version of your PMA.live
7. Edit and other options

The screenshot shows the 'Regions' page in the 'PMA.control test environment' v2.0.0.1338. The left sidebar contains navigation links: Home, Users, Projects, Slides & Images, Quizzes, Training sessions, Forms, System, Settings, Regions, System logs, Viewer settings, Speed test, Registration Terms, Attributes, Comment fields templates, Backend, and Documentation. The main content area is titled 'Regions' and includes a breadcrumb trail: Home / System settings / Regions. Below the title is a table with the following columns: Name, Geography, PMA.core URL, Information, PMA.live URL, and Information. A single instance is listed with the name 'PMA.core', location 'Strassburg', and version '2.1.0.2031'. The PMA.live URL is '2.0.0.98'. The page also includes a 'Create' button and an 'Edit' button.

1	2	3	4	5	6	7
Name	Geography	PMA.core URL	Information	PMA.live URL	Information	
PMA.core	Strassburg		2.1.0.2031 Valid license		2.0.0.98	Edit

You can also create a new instance by clicking the 'create' button to the top right of the table.

This will prompt you to fill in the Name, Location, and URLs of the PMA.core and PMA.live you wish to add.

Create new region

Name

Geography

PMA.core URL

PMA.live URL

Create

Cancel

By clicking 'Edit' (7) you will be given the same freetext fields to edit

Edit region

Home

System settings

Regions

Edit

EDIT REGION

Name

PMA.core

Geography

Strassburg

PMA.core URL

PMA.live URL

Save

Back to list

From:

<https://docs.pathomation.com/pathotrainer/> - **Pathotrainer**

Permanent link:

<https://docs.pathomation.com/pathotrainer/doku.php?id=regions>

Last update:

**2024/01/03 17:21**

