

# Regions

In the Regions section of the System section of the administrators toolbar, you will see your servers and geo-replication instances.

In this table you can see

1. The name of your instance
2. Where it is located
3. The web address of your PMA.core
4. The version and license information
5. The web address of your PMA.live
6. The version of your PMA.live
7. Edit and other options

The screenshot shows the 'Regions' section of the PMA control test environment. The interface includes a sidebar with navigation options like Home, Users, Projects, Slides & Images, Quizzes, Training sessions, Forms, System, Settings, and Backend. The main content area displays a table of regions with columns for Name, Geography, PMA.core URL, Information, PMA.live URL, and Information. A 'Create' button is visible in the top right of the table. The table contains one entry for 'PMA.core Strassburg' with version '2.1.0.2091' and license 'Valid license'. An 'Edit' button is also present for this entry.

You can also create a new instance by clicking the 'create' button to the top right of the table.

This will prompt you to fill in the Name, Location, and URLs of the PMA.core and PMA.live you wish to add.

### Create new region ✕

Name

Geography

PMA.core URL

PMA.live URL

By clicking 'Edit' (7) you will be given the same freetext fields to edit

### Edit region

[Home](#) / [System settings](#) / [Regions](#) / [Edit](#)

#### EDIT REGION

Name

Geography

PMA.core URL

PMA.live URL

From: <https://docs.pathomation.com/pathotrainer/> - **Pathotrainer**

Permanent link: <https://docs.pathomation.com/pathotrainer/doku.php?id=regions&rev=1704291541>

Last update: **2024/01/03 17:19**

