

What's new in version 2.1?

Pathotrainer 2.1.0 is a commit (1f204cca) of the Pathotrainer codebase (previously known as PMA.control). It follows [PMA.core 3.0.2](#) (be720dc8). It is a major revision, that contains a multiple new features and bug fixes.

New Features

Bugs fixed

- Fixed bug where the indexing mechanism could result in data loss.

Known issues

- Gamma parameter when passed by the client side, has no effect for CZI files(CORE3-68)
- Clicking on a user of a group returns an error (CORE-248)

See also

[PMA.core release history](#)

From:

<https://docs.pathomation.com/pathotrainer/> - **Pathotrainer**

Permanent link:

https://docs.pathomation.com/pathotrainer/doku.php?id=release_notes&rev=1715345834

Last update: **2024/05/10 15:57**

