

# What's new in version 2.1?

Pathotrainer 2.1.0 is a commit (1f204cca) of the Pathotrainer codebase (previously known as PMA.control). It follows \_\_. It is a major release that contains a multiple new features and bug fixes.

## New Features

- CTRL-96: Users can now make annotations in a session. These are only visible to supervisors, and the supervisor is able to filter them
- CTRL-97: Sharing of WSIs and ROIs is now possible
- CTRL-98: Project owners / Supervisors can now set how all image labels will be displayed ( i.e. filenames | captions | no label ) with a single click in the Training Session options.

## Bugs fixed

- CTRL-145: Fixed bug in which WSI content would not load in Course Mode

## Known issues

- Gamma parameter when passed by the client side, has no effect for CZI files(CORE3-68)
- Clicking on a user of a group returns an error (CORE-248)

## See also

[PMA.core release history](#)

From:  
<https://docs.pathomation.com/pathotrainer/> - **Pathotrainer**

Permanent link:  
[https://docs.pathomation.com/pathotrainer/doku.php?id=release\\_notes&rev=1715346803](https://docs.pathomation.com/pathotrainer/doku.php?id=release_notes&rev=1715346803)

Last update: **2024/05/10 16:13**

