

Revealing Annotations

Prior to August 2025 it was possible to give immediate feedback to independent learners completing cases in the Self-Test interaction mode in the form of the correct diagnosis and additional text comments.

With the release of PathoTrainer 2.1.0, it is now possible to also guide your students with on-slide drawn annotations that reveal themselves upon score submission

See how this works, and learn how to apply these yourself in the Consensus Score Edit mode in the instructional video below:

[pathotrainer_revealing_annotations.mp4](#)

From:

<https://docs.pathomation.com/pathotrainer/> - **Pathotrainer**

Permanent link:

https://docs.pathomation.com/pathotrainer/doku.php?id=reveal_expert_annotations

Last update: **2025/11/26 18:05**

