

Session page

Once you join a session, the system takes you to the session menu. On this page you will find all material related to the session. The exact layout of a session may vary depending on the configuration, with more or less options in each section.

The sections include:

1. Introductory material - this can include instructions/user guides, or questionnaires to fill in before completing course content. If you have completed a form, a greentext 'Done' will be displayed.
2. Case collections - these can appear in different states e.g. locked, active test, browsing mode etc.. depending on how the supervisors have currently configured the training session. If a collection is available to score - as is the case for Test 3 in the example below - the number of cases you have evaluated will be displayed in blue text
3. Your overall score
4. Supplemental material e.g. videos and website links. The screen compatibility check will always be present in every session, and allows you to quickly check if the contrast of your screen is sufficient for digital pathology.
5. Contact information of supervisors and administrators who you can reach out to if you have questions or problem.

PMA.control test environment

Training Session Test

- 1 Before starting, please take a few moments to complete a short questionnaire
? Test form Done ✓
- 2 Case collections
A project to show how easy it is to start a project!
Test Compare the consensus score to yours
Test 3 0/9 cases evaluated Take the test by scoring each case
Test 2 Compare the consensus score to yours
- 3 My score
- 4 Supplemental material
Screen suitability test
- 5 Do you need help?
Sea Bow - @pathomation.com
Project owners
Sea Bow - @pathomation.com

From:

<https://docs.pathomation.com/pathotrainer/> - **Pathotrainer**

Permanent link:

https://docs.pathomation.com/pathotrainer/doku.php?id=session_main_page

Last update: **2024/01/03 15:57**

