


## Session main page (supervisor)

Once you join a session, the system takes you to the main page of it. Here you can access the available case collections as well as supplemental material related to the session. The exact layout of a session may vary depending on the configuration. A typical setup is illustrated below.



In addition to the material presented to the trainees/participants a trainer can also get access to additional material that isn't presented to the trainees. E.g. documents containing consensus-scores - additional comments etc.

A trainer also has access to:

- the participants management page
- comprehensive score report of all participants.

From:  
<https://docs.pathomation.com/pathotrainer/> - **Pathotrainer**

Permanent link:  
[https://docs.pathomation.com/pathotrainer/doku.php?id=session\\_main\\_page\\_supervisor&rev=1597757743](https://docs.pathomation.com/pathotrainer/doku.php?id=session_main_page_supervisor&rev=1597757743)

Last update: **2020/08/18 16:35**

