

Setting Up a Session

Previously if you wanted to work with slides for a biomarker training or consensus scoring you'd need a full briefcase with exemplar slides for all your participants per indication/biomarker, forms, case information... and that's before you even begin to consider the microscopes!

Training Sessions were designed to completely eradicate these difficult and costly logistical hurdles. A standard training session contains everything you need to train other pathologists, reach a ground-truth consensus, collect annotations from experts around the globe, or validate an AI algorithm - ask pm@pathomation.com for information on any of these or to outsource setting these up.



A Training Session contains:

1. All your participants
2. The ability to invite and/or contact all of them at once
3. Tracking of email invitations and logins
4. Forms to capture the data from the participants
5. Supplementary info e.g. manuals / interpretation guides (videos, powerpoints, pdf's, etc)
6. Case collections ('slide boxes') which all your participants can see at once
7. Consensus scores
8. Optional Annotations made either by the Supervisor, and/or by the Trainee Participants

Walkthrough

This video will get you started:

[pathotrainer_setup_training_session.mp4](#)

You will also see a bluetext prompt to create a Training Session when viewing a Project that does not yet have a Training Session

The screenshot shows the Pathotrainer web application. On the left is a dark sidebar with a menu: Home, Users, Projects (highlighted with a red box), Slides & Images, Quizzes, Training sessions, Forms, System, Backend, and Documentation. The main area displays the 'Patho_Project' details. At the top, the breadcrumb 'Home / Projects / Patho_Project' is visible, with 'Patho_Project' highlighted in a red box. Below the breadcrumb are tabs: Details (selected), Edit, Additional material, Links, Integrity check, and Annotations toolbar. The project details include: Title (Patho_Project), Logo, State (Active), Description (A project to show how easy it is to start a project!), Contact info, Viewer theme (Modern), and Project Owners. A section for 'Case collections [1]' contains a table:

Title	# of cases	Scoring Form
Test	1	Test form

Below the table, under 'Training sessions', it says 'There are no training sessions in this project.' and includes a link 'Setup a training session' which is pointed to by a red arrow.

Well done! You've created your first training session. Now it's time to [add participants](#) or [put on some finishing touches](#)

From:

<https://docs.pathomation.com/pathotrainer/> - **Pathotrainer**

Permanent link:

https://docs.pathomation.com/pathotrainer/doku.php?id=setup_a_session

Last update: **2025/11/20 13:40**

