

As a trainer-supervisor you are a super-user in a session. You will be able to change the way the participants will see the slides. Depending on the predefined interaction modes participants will be able to see slides with empty dataforms or complete score overviews. As a trainer-supervisor you can also take over all the screens of the participants and start an interactive live session.

From:

<https://docs.pathomation.com/pathotrainer/> - **Pathotrainer**

Permanent link:

<https://docs.pathomation.com/pathotrainer/doku.php?id=supervisor-trainer&rev=1635363957>

Last update: **2021/10/27 22:45**

