

Let's break down the components to setting up a session. They are summarized below.

- select the case collections and data forms that will be used
- edit the session welcome page all the participants will see.
- decide whether you will send email invitations
- add all participants to the session and send them invites from pma.control software.

From:

<https://docs.pathomation.com/pathotrainer/> - **Pathotrainer**

Permanent link:

[https://docs.pathomation.com/pathotrainer/doku.php?id=what\\_is\\_a\\_session&rev=1635701912](https://docs.pathomation.com/pathotrainer/doku.php?id=what_is_a_session&rev=1635701912)

Last update: **2021/10/31 20:38**

