


Conference settings (administration)

This page allows you to setup PMA.Live server URL that creates live broadcast session for a slide to be viewed by audience using PMA.Studio. You can also validate the PMA.Live URL before saving.



The screenshot shows the Pathomation administration interface. On the left is a dark sidebar with a menu containing: Title Servers, Security, Panels Layout, Ribbon Layout, Global Settings, Shared Links, Conference Settings (highlighted in green), Server Visits, Users, Manual, Swagger Docs, and Viewer. The main content area is titled 'Configure Conferencing'. It features a 'PMA.Live uri' label, a text input field containing 'https://host.pathomation.com/pma/live/', a 'Validate' button, and 'Save' and 'Reset' buttons below the input field.

From:

<https://docs.pathomation.com/pimsls/3.0.0/> - **PIMS LS 3.0.0**

Permanent link:

https://docs.pathomation.com/pimsls/3.0.0/doku.php?id=admin_conference

Last update: **2023/04/06 22:58**

