

## Shared links (administration)

This page list down all the shared links created for a slide or slide folder. It provides detail about each of the shared links that include link type, is password protected and expiry date. You can also delete or suspend a shared link.

| Shared Links |             |         |       |                          |                                        |                                                                                     |                                                                                     |
|--------------|-------------|---------|-------|--------------------------|----------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| Id           | Name        | Comment | Type  | Is protected by password | Expires At                             | Actions                                                                             |                                                                                     |
| 03U6dma      |             |         | ROI   | Not Protected            | Jan 01 0001 00:00 <span>Expired</span> |  |  |
| 04qVAKUj     | testje Cleg |         | Slide | Not Protected            | Never                                  |  |  |
| 05mKU1UQ     |             |         | Slide | Not Protected            | Never                                  |  |  |
| 0UGVFUj      |             |         | Grid  | Not Protected            | Never                                  |  |  |
| 0JlmmT       |             |         | Slide | Not Protected            | Never                                  |  |  |
| 0SHKCT7aY    | Cleg test   |         | ROI   | Protected                | Jul 08 2021 21:00 <span>Expired</span> |  |  |
| 0w26CH9      |             |         | ROI   | Not Protected            | Never                                  |  |  |
| 1Aqpsdq      |             |         | ROI   | Not Protected            | Never                                  |  |  |
| 1136dHMQ     |             |         | Slide | Not Protected            | Aug 17 2020 19:00 <span>Expired</span> |  |  |
| 1NC933W      |             |         | Slide | Not Protected            | Never                                  |  |  |
| 1KPmwS1E     |             |         | Grid  | Not Protected            | Never                                  |  |  |
| 1M3F0My      |             |         | Slide | Not Protected            | Never                                  |  |  |

From:

<https://docs.pathomation.com/pimsls/3.0.0/> - **PIMS LS 3.0.0**

Permanent link:

[https://docs.pathomation.com/pimsls/3.0.0/doku.php?id=admin\\_links](https://docs.pathomation.com/pimsls/3.0.0/doku.php?id=admin_links)

Last update: **2023/04/06 22:57**

