

## Available ribbon commands

| Ribbon section | Command name               | Parameters |
|----------------|----------------------------|------------|
| File menu      | custom-panel-algo-ichi     | props=     |
|                | custom-panel-algo-ni       | props=     |
|                | file                       |            |
|                | about-us                   |            |
|                | show-eula                  |            |
|                | manual                     |            |
|                | sign-out                   |            |
| Home tab       | simple-layout              |            |
|                | detailed-layout            |            |
|                | detailed-layout            |            |
|                | viewport-toggle-overview   |            |
|                | viewport-toggle-slide-path |            |
|                | viewport-toggle-barcode    |            |
|                | add-tab                    |            |
|                | reset-layout               |            |
|                | viewport-toggle-zoom       |            |
|                | viewport-toggle-scalebar   |            |
|                | viewport-toggle-label      |            |
|                | panel-toggle-color         |            |
|                | panel-toggle-search        |            |
|                | open-folder-in-grid        |            |
|                | prev-directory             |            |
|                | next-directory             |            |
|                | prev-image                 |            |
|                | next-image                 |            |
|                | add-folder-bookmark        |            |
|                | reset-grid                 |            |
|                | clear-cell                 |            |
|                | open-in-tab                |            |
|                | popout-info                |            |
|                | snapshot                   |            |
|                | quick-snapshot             |            |
|                | toggle-toolbar             |            |
|                | increase-zoom              |            |
|                | decrease-zoom              |            |
|                | custom-zoom-1x             |            |
|                | custom-zoom-2x             |            |
|                | custom-zoom-5x             |            |
|                | panel-toggle-color         |            |
|                | send-slide-to-tray         |            |
|                | send-to-tray               |            |
|                | panel-toggle-tray          |            |
|                | share-slide-tray           |            |

| Ribbon section  | Command name                                | Parameters |
|-----------------|---------------------------------------------|------------|
|                 | <a href="#">quick-share-folder</a>          |            |
|                 | <a href="#">quick-share-slide</a>           |            |
|                 | <a href="#">quick-share-slide-roi</a>       |            |
|                 | <a href="#">quick-share-grid</a>            |            |
|                 | <a href="#">quick-share-grid-rois</a>       |            |
|                 | <a href="#">quick-share-thumbnail</a>       |            |
|                 | <a href="#">share-slide</a>                 |            |
| Grid tab        | <a href="#">toggle-1x1-grid</a>             |            |
|                 | <a href="#">toggle-2x2-grid</a>             |            |
|                 | <a href="#">toggle-3x3-grid</a>             |            |
|                 | <a href="#">sync-views-toggle</a>           |            |
|                 | <a href="#">reset-grid</a>                  |            |
|                 | <a href="#">clear-cell</a>                  |            |
|                 | <a href="#">open-in-tab</a>                 |            |
|                 | <a href="#">viewport-add-row</a>            |            |
|                 | <a href="#">viewport-remove-row</a>         |            |
|                 | <a href="#">viewport-remove-row</a>         |            |
|                 | <a href="#">viewport-remove-column</a>      |            |
|                 | <a href="#">flip-rotate-menu</a>            |            |
|                 | <a href="#">reset-flip</a>                  |            |
|                 | <a href="#">flip-horizontally</a>           |            |
|                 | <a href="#">flip-vertically</a>             |            |
|                 | <a href="#">reset-rotation</a>              |            |
|                 | <a href="#">rotate-clockwise-90</a>         |            |
|                 | <a href="#">rotate-counter-clockwise-90</a> |            |
|                 | <a href="#">rotate-180</a>                  |            |
|                 | <a href="#">zoom-in</a>                     |            |
|                 | <a href="#">panel-toggle-annotations</a>    |            |
|                 | <a href="#">save-annotations</a>            |            |
|                 | <a href="#">toggle-annotations</a>          |            |
|                 | <a href="#">toggle-annotation-text</a>      |            |
|                 | <a href="#">toggle-fingerprint</a>          |            |
| Annotations tab | <a href="#">Cursor</a>                      |            |
|                 | <a href="#">delete-annotation</a>           |            |
|                 | <a href="#">delete-all-annotations</a>      |            |
|                 | <a href="#">edit-annotation-text</a>        |            |
|                 | <a href="#">merge-selected-annotations</a>  |            |
|                 | <a href="#">modify-selected-annotation</a>  |            |
|                 | <a href="#">Point</a>                       |            |
|                 | <a href="#">Circle</a>                      |            |
|                 | <a href="#">Line</a>                        |            |
|                 | <a href="#">Arrow</a>                       |            |
|                 | <a href="#">Ellipse</a>                     |            |
|                 | <a href="#">Rectangle</a>                   |            |
|                 | <a href="#">Landmark1</a>                   |            |

| Ribbon section   | Command name                       | Parameters                               |
|------------------|------------------------------------|------------------------------------------|
|                  | MultiPoint                         |                                          |
|                  | Polygon                            |                                          |
|                  | ClosedFreehand                     |                                          |
|                  | CompoundFreehand                   |                                          |
|                  | Freehand                           |                                          |
|                  | Wand                               |                                          |
|                  | Brush                              |                                          |
|                  | export-to-wkt                      |                                          |
|                  | export-to-mld                      |                                          |
|                  | export-to-annotations              |                                          |
|                  | export-to-aperio                   |                                          |
|                  | export-to-csv                      |                                          |
|                  | import-to-wkt                      |                                          |
|                  | import-to-mld                      |                                          |
|                  | import-to-halo                     |                                          |
|                  | import-to-csv                      |                                          |
| Pre-set tab      | preset-Circle0237mm2               | pma-shape, pma-size-a, pma-size-b        |
|                  | preset-Circle1mm                   | pma-shape, pma-size-a, pma-size-b        |
|                  | preset-Circle2mm                   | pma-shape, pma-size-a, pma-size-b        |
|                  | preset-Circle1cm                   | pma-shape, pma-size-a, pma-size-b        |
|                  | preset-Square5mm                   | pma-shape, pma-size-a, pma-size-b        |
|                  | preset-Square1cm                   | pma-shape, pma-size-a, pma-size-b        |
|                  | preset-Tumor                       | pma-classification, pma-shape, pma-color |
|                  | preset-Macro                       | pma-classification, pma-shape, pma-color |
| Conferencing tab | create-conference                  |                                          |
|                  | join-conference                    |                                          |
|                  | start-conference                   |                                          |
|                  | pause-conference                   |                                          |
|                  | stop-conference                    |                                          |
|                  | toggle-everyone-control-conference |                                          |
|                  | broadcast-conference               |                                          |
|                  | panel-toggle-conference            |                                          |
|                  | panel-toggle-forms                 |                                          |
|                  | panel-toggle-links                 |                                          |
|                  | panel-toggle-navigation            |                                          |
|                  | panel-toggle-bookmarks             |                                          |

From:

<https://docs.pathomation.com/pimsls/3.0.0/> - PIMS LS 3.0.0

Permanent link:

[https://docs.pathomation.com/pimsls/3.0.0/doku.php?id=ribbon\\_commands&rev=1674030624](https://docs.pathomation.com/pimsls/3.0.0/doku.php?id=ribbon_commands&rev=1674030624)

Last update: 2023/01/18 11:30

