

Ribbon > Home tab > Layout group > Custom layout

Button identification



Associated command

panels-dialog

Function description

Brings up a dialog that offers the end-user fine-grained control over what [panels](#) are visible:

Items	Show
Slide browser	☒
Forms Data	☐
Slide Info	☒
Annotations	☒
Tray	☐
Conference	☐
Color Adjustment	☐
Shared Links	☐
Search panel	☐
Bookmarks panel	☐
Add grid	☒
04310688C0010S.scn	☒
Add IFrame	☒

Related functionality

[Ribbon > Home tab](#)

From:

<https://docs.pathomation.com/pimsls/3.0.0/> - PIMS LS 3.0.0

Permanent link:

https://docs.pathomation.com/pimsls/3.0.0/doku.php?id=ribbon_home_layout_custom

Last update:

2023/03/03 16:25



