

Ribbon > Home tab > Layout group > Default layout

Button identification




Associated command

simple-layout

Function description

This restores the panel layout to a two-panel layout with the [slide explorer panel](#) on a left and a [viewport panel](#) on the right.



This button always resets to the same layout, irrespective of the application's default layout as stipulated in the [Administrative interface to define default panel configuration](#).

Related functionality

[Ribbon > Home tab](#)

From:

<https://docs.pathomation.com/pimsls/3.0.0/> - **PIMS LS 3.0.0**

Permanent link:

https://docs.pathomation.com/pimsls/3.0.0/doku.php?id=ribbon_home_layout_default

Last update: **2023/02/10 18:41**

