

Ribbon > Home tab > Layout group > Default layout


Button identification



Associated command

Function description

This restores the panel layout to a two-panel layout with the [slide explorer panel](#) on the left and a [viewport panel](#) on the right.



Related functionality

[Ribbon > Home tab](#)

Last
update: 2023/02/08 15:44 ribbon_home_layout_default https://docs.pathomation.com/pimsls/3.0.0/doku.php?id=ribbon_home_layout_default&rev=1675860240

From:
<https://docs.pathomation.com/pimsls/3.0.0/> - **PIMS LS 3.0.0**

Permanent link:
https://docs.pathomation.com/pimsls/3.0.0/doku.php?id=ribbon_home_layout_default&rev=1675860240

Last update: **2023/02/08 15:44**

