


Grid: Slide Tool Group

The Slide dropdown contains tools for turning changing brightness/contrast of your current image, rotating and magnifying the slide.



Clear current cell - clears current cell of the grid Clear all grid cells - clears all cells of a grid Open slide from current cell into a new viewport panel Popout selected slide in a new window

From:

<https://docs.pathomation.com/pimsls/3.0.0/> - **PIMS LS 3.0.0**

Permanent link:

<https://docs.pathomation.com/pimsls/3.0.0/doku.php?id=slide&rev=1746905182>

Last update: **2025/05/10 22:26**

