

Interacting with slides


PMA.studio provides multiple ways to interact with slides that includes:

1. Working with grids
2. Static snapshots
3. Sharing live views
4. Broadcasting and conferencing

UI provides tabs “GRID” and “CONFERENCING” that allows you to interact with slides and perform different actions on it.

Grid Tab


Grid tab provides set of tools that allows resizing grid size, slide static snapshot and live sharing of slides.



[See more details on what each button on the grid tab does.](#)

(Conferencing Tab)

This unique tab allows you to create a conference or to broadcast a live session where an expert can add annotations to a slide to be viewed by number of audience.



From:

<https://docs.pathomation.com/pimsls/3.0.0/> - **PIMS LS 3.0.0**

Permanent link:

<https://docs.pathomation.com/pimsls/3.0.0/doku.php?id=slides&rev=1677854695>

Last update: **2023/03/03 17:44**

