


Working with grids

PMA.studio uses viewport as a grid system. UI provides a separate tab “GRID” to perform grid operations on the viewport.




“GRID” tab consists of 6 sections:

1. Size
2. Synchronize
3. R/C
4. Slide
5. Snapshot
6. Share


Size

Size section allows you to break viewport as a grid system with different sizes. This includes grid sizes of 1×1, 2×2, 3×3 and even you can customized the grid size according to your need.




(1×1 Grid size)

First grid option allows you to divide the viewport into 1×1 column grid.




(2×2 Grid size)

Second grid option divides the viewport into 2×2 column grid.




(3x3 Grid size)

Third grid option divides the viewport into 3x3 column grid.




(Custom Grid size)

You can even customize the grid according to your need by using the custom grid option.



Synchronize

This section allows you to perform different actions as a whole on individual grid column or as a whole grid system. It is further divided into 4 sub actions:



1. Toggle synchronous navigation between cells in the current grid
2. * Unordered List Item
3. Clear all grid cells
4. Clear current cell only
5. Open current slide into separate viewport panel

From:

<https://docs.pathomation.com/pimsls/3.0.0/> - **PIMS LS 3.0.0**

Permanent link:

https://docs.pathomation.com/pimsls/3.0.0/doku.php?id=slides_grid&rev=1650996093

Last update: **2022/04/26 21:01**

