

ELMAH (Error Logging Modules and Handlers) is an application-wide error logging facility that is completely pluggable. It captures nearly all unhandled exceptions in PMA.studio as well as some handled exceptions that needs to be reported. It also provides the following very useful features for maintaining a running production server:

- Logging of nearly all unhandled exceptions.
- A web page to remotely view the entire log of recorded exceptions.
- A web page to remotely view the full details of any one logged exception, including colored stack traces.

Elmah although powerful and useful it is, it is not a substitute for IIS logs and PMA.studio internal logging as it does not capture all detailed information for requests like IIS does, nor does know the internal workings of PMA.studio as the internal logging is capable of. Once you click on the “Elmah Logs” link in admin panel you are taken to the errors page. You can click on the “detail” link to view more details about the specific error.

ELMAHCore

RSS FeedRSS ETagDownload LogHelpAbout

Error Log for PMA.Studio on MALTA1767

Errors 1 to 15 of total 96 (page 1 of 7). Start with [15](#), [16](#), [20](#), [25](#), [30](#), [50](#) or [100](#) errors per page.

Host	Code	Type	Server	User	When
MALTA1767	404	HTTP	NotFound Details...		yesterday
MALTA1767	404	HTTP	NotFound Details...		yesterday
MALTA1767	404	HTTP	NotFound Details...		yesterday
MALTA1767	404	HTTP	NotFound Details...		yesterday
MALTA1767	404	HTTP	NotFound Details...		yesterday
MALTA1767	404	HTTP	NotFound Details...		yesterday
MALTA1767	404	HTTP	NotFound Details...		yesterday
MALTA1767	404	HTTP	NotFound Details...		yesterday
MALTA1767	404	HTTP	NotFound Details...		yesterday
MALTA1767	404	HTTP	NotFound Details...		yesterday
MALTA1767	404	HTTP	NotFound Details...		yesterday
MALTA1767	404	HTTP	NotFound Details...		yesterday
MALTA1767	404	HTTP	NotFound Details...		yesterday
MALTA1767	404	HTTP	NotFound Details...		yesterday
MALTA1767	404	HTTP	NotFound Details...		yesterday

Next errors

Based on [Fluent](#). Downloaded by [FluentCase](#), version 1.2.6.8. [FluentCase](#) is licensed under [Apache License, Version 2.0](#). Server date is Thursday, 04 August 2022. Server time is 17:51:17. All dates and times displayed are in the UTC+0 (GMT) Time zone. This box is provided by the Microsoft SQL Server Profiler | [View all logs...](#)

Based on [Elastic](#). Powered by [Dinohq/elastic](#), version 1.2.6.8. Dinohq/elastic Licensed under [Apache License, Version 2.0](#). Server date is Thursday, 04 August 2023. Server time is 17:31:17. All dates and times displayed are in the W. Europe Daylight Time zone. This too is provided by the [Microsoft SQL Server Error Log](#).

Last update: **2022/08/04 18:37**

