

Global settings (administration)


This admin page allows us to set global settings for PMA.Studio application. Page is divided into two set of fields:

Global settings has 5 global fields: ·

1. Language
2. Application IconNo file chosen
3. Email For Help
4. Conference Server Url
5. Snapshot Filename Template

Annotation settings has 5 global fields: ·

1. Show Annotation Dialog
2. Auto Save Annotations
3. Restart Annotations Drawing
4. Show Slide Info Button
5. Show Full Slide Path In Viewport



From:

<https://docs.pathomation.com/pimsls/3.0.1/> - PIMS LS 3.0.1

Permanent link:

https://docs.pathomation.com/pimsls/3.0.1/doku.php?id=admin_global&rev=1648660454

Last update: **2022/03/30 20:14**

