

Global settings (administration)


This admin page allows us to set global settings for PMA.Studio application. Page is divided into two set of fields:

Global settings has 6 fields:

1. Language
2. Application IconNo file chosen
3. Email For Help
4. Conference Server Url
5. Snapshot Filename Template

Annotation settings has 5 fields:

1. Show Annotation Dialog
2. Auto Save Annotations
3. Restart Annotations Drawing
4. Show Slide Info Button
5. Show Full Slide Path In Viewport



From:
<https://docs.pathomation.com/pimsls/3.0.1/> - **PIMS LS 3.0.1**

Permanent link:
https://docs.pathomation.com/pimsls/3.0.1/doku.php?id=admin_global&rev=1680811111

Last update: **2023/04/06 22:58**

