

===== Tile servers (administration) =====c

In this section of PMA.studio's administrative portal, you can manage the different PMA.core instances that are used by PMA.studio to retrieve content from.

PMA.studio (PMA.view) does not function without connectivity to a PMA.core tile server. If you want to learn more about what a tile server is, [we have a YouTube video on this topic](#).

From:

<https://docs.pathomation.com/pimsls/3.0.1/> - **PIMS LS 3.0.1**

Permanent link:

[https://docs.pathomation.com/pimsls/3.0.1/doku.php?id=admin\\_tileservers&rev=1643046669](https://docs.pathomation.com/pimsls/3.0.1/doku.php?id=admin_tileservers&rev=1643046669)

Last update: **2022/01/24 20:51**

