


Server visits (administration)

This page list down all the different visits to PMA.Studio. Records can be filtered through PMA.core server url, username and date of visit. Each row consists of:

1. PMA.core server url accessed
2. PMA.Studio username used for login
3. Date of visit
4. Status of success/failure on login
5. IP address through which the user connected
6. Browser/Device used to login



Server URL	Username	Date	Status	IP Address	User Agent
https://docs.pathomation.com/pimsls/3.0.1/	h_j_pert	Wednesday, March 29th 2023, 11:45:53 am	Success	208.84.186.75	Mozilla/5.0 (Windows NT 10.0; Win64; x64; AppleWebKit/537.36; Chrome/109.0.5414.106) Safari/537.36
https://docs.pathomation.com/pimsls/3.0.1/	h_j_pert	Wednesday, March 29th 2023, 12:03:16 pm	Success	208.84.186.75	Mozilla/5.0 (Windows NT 10.0; Win64; x64; AppleWebKit/537.36; Chrome/109.0.5414.106) Safari/537.36
https://docs.pathomation.com/pimsls/3.0.1/	h_j_pert	Wednesday, March 29th 2023, 12:03:40 pm	Success	208.84.186.75	Mozilla/5.0 (Windows NT 10.0; Win64; x64; AppleWebKit/537.36; Chrome/109.0.5414.106) Safari/537.36
https://docs.pathomation.com/pimsls/3.0.1/	h_j_pert	Wednesday, March 29th 2023, 11:53:54 am	Success	208.84.186.75	Mozilla/5.0 (Windows NT 10.0; Win64; x64; AppleWebKit/537.36; Chrome/109.0.5414.106) Safari/537.36
https://docs.pathomation.com/pimsls/3.0.1/	h_j_pert	Wednesday, March 29th 2023, 12:45:40 pm	Success	208.84.186.75	Mozilla/5.0 (Windows NT 10.0; Win64; x64; AppleWebKit/537.36; Chrome/109.0.5414.106) Safari/537.36
https://docs.pathomation.com/pimsls/3.0.1/	h_j_pert	Thursday, March 29th 2023, 12:20:40 pm	Success	208.84.186.75	Mozilla/5.0 (Windows NT 10.0; Win64; x64; AppleWebKit/537.36; Chrome/109.0.5414.106) Safari/537.36
https://docs.pathomation.com/pimsls/3.0.1/	h_j_pert	Thursday, March 29th 2023, 12:17:17 pm	Success	208.84.186.75	Mozilla/5.0 (Windows NT 10.0; Win64; x64; AppleWebKit/537.36; Chrome/109.0.5414.106) Safari/537.36
https://docs.pathomation.com/pimsls/3.0.1/	h_j_pert	Thursday, March 29th 2023, 12:18:33 pm	Success	208.84.186.75	Mozilla/5.0 (Windows NT 10.0; Win64; x64; AppleWebKit/537.36; Chrome/109.0.5414.106) Safari/537.36

From:

<https://docs.pathomation.com/pimsls/3.0.1/> - PIMS LS 3.0.1

Permanent link:

https://docs.pathomation.com/pimsls/3.0.1/doku.php?id=admin_visits&rev=1648666482

Last update: **2022/03/30 21:54**

