


# Annotations

PIMS LS is the swiss army knife of WSI annotations - besides the extensive annotation tools provided, it also supports many formats of annotations made in other programs.



So, that was just the beginning and we have not many any actual annotation as yet. Now, if we look at the “Annotation” tab content on the ribbon it consists of different types of annotations tools we can use to draw annotations on a selected slide.

PMA.Studio provides different set of annotations tools you can choose from and draw annotations:

## (Simple shapes)

1. Point
2. Circle
3. Line
4. Arrow
5. Ellipse
6. Rectangle
7. Text
8. Landmark 1
9. Landmark 2
10. Landmark 3

## (Polygon Shapes)

1. Multipoint
2. Polygon
3. Closed Freehand
4. Compound Freehand
5. Freehand
6. Wand
7. Brush

From:  
<https://docs.pathomation.com/pimsls/3.0.1/> - **PIMS LS 3.0.1**

Permanent link:  
<https://docs.pathomation.com/pimsls/3.0.1/doku.php?id=annotations&rev=1732276296>

Last update: **2024/11/22 14:51**

