

# Brush and magic wand


Brush and magic wand are two special tools that allow you to draw annotation on specific area of a slide.

(Brush Tool)

Brush tool allows you to paint specific area and create annotation on it.

(Wand Tool)

Wand tool allows you to select different areas on a slide by simply clicking on the slide. Once you click on a slide that specific area is automatically selected and annotation can be created on that area by clicking “create annotation from selection” button.



From:

<https://docs.pathomation.com/pimsls/3.0.1/> - PIMS LS 3.0.1

Permanent link:

[https://docs.pathomation.com/pimsls/3.0.1/doku.php?id=annotations\\_magic&rev=1650563012](https://docs.pathomation.com/pimsls/3.0.1/doku.php?id=annotations_magic&rev=1650563012)

Last update: 2022/04/21 20:43

