


# Configure Viewport

The Configure Viewport workgroup contains icons for toggling on/off different viewport elements



In the picture above you can see the 6 toggles and the parts of the UI they control:

1. Overview
2. File path
3. Barcode (this one is anonymized)
4. Zoom control
5. Scale bar
6. Slide label

From:

<https://docs.pathomation.com/pimsls/3.0.1/> - **PIMS LS 3.0.1**

Permanent link:

<https://docs.pathomation.com/pimsls/3.0.1/doku.php?id=configviewport>

Last update: **2024/12/04 17:18**

