


## Forms panel

The “Forms data” panel purpose is to show filled forms for the slide and provide access to editing and filling new forms for currently selected slide.

The panel list downs existing forms and also provide option to add/edit forms. To add a form user can click on the “Add” link on the right hand corner which will popup dialog box that shows up list of forms to choose from. Once user selects a form and press “add” button it loads the respective form to the panel.



From:

<https://docs.pathomation.com/pimsls/3.0.1/> - PIMS LS 3.0.1

Permanent link:

[https://docs.pathomation.com/pimsls/3.0.1/doku.php?id=panel\\_forms&rev=1680722577](https://docs.pathomation.com/pimsls/3.0.1/doku.php?id=panel_forms&rev=1680722577)

Last update: **2023/04/05 22:22**

