

The Forms Panel

The Forms panel is the place to input and review additional data about your slides.

You can view a quick video guide on how to use the forms panel below:

[pimsls_formshowto.mp4](#)

Please note forms need to be defined / created in PMA.core. Please ask your administrator if you don't have the form(s) you need.

From:
<https://docs.pathomation.com/pimsls/3.0.1/> - **PIMS LS 3.0.1**

Permanent link:
https://docs.pathomation.com/pimsls/3.0.1/doku.php?id=panel_forms&rev=1733481448

Last update: **2024/12/06 13:37**

