

## Ribbon > Home tab > Layout group > Custom layout

### Button identification



### Associated command

panels-dialog

### Function description

Brings up a dialog that offers the end-user fine-grained control over what [panels](#) are visible:

✚

Items	Show
Slide browser	<input checked="" type="checkbox"/>
Forms Data	<input type="checkbox"/>
Slide Info	<input checked="" type="checkbox"/>
Annotations	<input checked="" type="checkbox"/>
Tray	<input type="checkbox"/>
Conference	<input type="checkbox"/>
Color Adjustment	<input type="checkbox"/>
Shared Links	<input type="checkbox"/>
Search panel	<input type="checkbox"/>
Bookmarks panel	<input type="checkbox"/>
Add grid	<input checked="" type="checkbox"/>
04310688C0010S.scn	<input checked="" type="checkbox"/>
Add IFrame	<input checked="" type="checkbox"/>

Close

## Related functionality

[Ribbon > Home tab](#)

From: <https://docs.pathomation.com/pimsls/3.0.1/> - PIMS LS 3.0.1

Permanent link: [https://docs.pathomation.com/pimsls/3.0.1/doku.php?id=ribbon\\_home\\_layout\\_custom](https://docs.pathomation.com/pimsls/3.0.1/doku.php?id=ribbon_home_layout_custom)

Last update: 2023/03/03 16:25



