

Ribbon > Home tab > Layout group > Default layout

Button identification




Associated command

simple-layout

Function description

This restores the panel layout to a two-panel layout with the [slide explorer panel](#) on a left and a [viewport panel](#) on the right.



Related functionality

[Ribbon > Home tab](#)

Last
update: 2023/02/10 16:31 ribbon_home_layout_default https://docs.pathomation.com/pimsls/3.0.1/doku.php?id=ribbon_home_layout_default&rev=1676035914

From:
<https://docs.pathomation.com/pimsls/3.0.1/> - **PIMS LS 3.0.1**

Permanent link:
https://docs.pathomation.com/pimsls/3.0.1/doku.php?id=ribbon_home_layout_default&rev=1676035914

Last update: **2023/02/10 16:31**

