



## Rotation Tool

Rotation dropdown consists of different rotation options for the viewport including flipping viewport horizontally, vertically, rotating it clockwise or counter clockwise and then reset the flips and rotations.



- Flip horizontally - flips viewport slide horizontally
- Flip vertically - flips viewport slide vertically
- Rotate counter-clockwise 90 degree
- Rotate 180 degree
- Reset flip
- Reset rotation



From:

<https://docs.pathomation.com/pimsls/3.0.1/> - PIMS LS 3.0.1

Permanent link:

<https://docs.pathomation.com/pimsls/3.0.1/doku.php?id=rotation>Last update: **2025/05/13 11:21**