



## Rotation Tools

Rotation dropdown consists of different rotation options for the viewport including flipping viewport horizontally, vertically, rotating it clockwise or counter clockwise and then reset the flips and rotations.



- Flip horizontally - flips viewport slide horizontally
- Flip vertically - flips viewport slide vertically
- Rotate counter-clockwise 90 degree
- Rotate 180 degree
- Reset flip
- Reset rotation



From:  
<https://docs.pathomation.com/pimsls/3.0.1/> - PIMS LS 3.0.1

Permanent link:  
<https://docs.pathomation.com/pimsls/3.0.1/doku.php?id=rotation&rev=1747059890>

Last update: 2025/05/12 17:24

