

Terminology and glossary


As with any technology, there are lots of terms we toss around to explain things. Just understanding what each term means can be half the battle. Whether you are entirely new to PMA.studio and PMA.vue or have been using it for years, there are always new terms to learn. The fact that we use common English words for many capabilities can add an additional layer of confusion.

This page is an attempt to demystify some of the terms and acronyms we use every day in working with the platform. See one missing? Feel free to [Contact us](#) to tell us what you would want to see added.

Area

Area is a group of pages inside web-application, PMA.Studio has 2 areas - Admin and Viewer. Admin area contains pages that control application state and configuration and viewer area consists of pages that user spent most of the time that gives access to slide viewing and related functionality.

Area for end-users



Area for administrators



Panels and panel stacks

A panel is a separate application component of specific type. There are n panels - Viewer, Tree, Tray and others, you can find more details in Panels section.

A panel stack is a group of panels organized in tabbed view each tab contains single panel and at one time only single tab can be visible in panel stack.

Ribbon

Ribbon is a group of controls (buttons, dropdowns etc.) placed on top of viewer area. On the illustration below you can see a Ribbon with specific configuration, it can be different on your side.

Tile server

A tile server is a web application - [PMA.core](#) or [PMA.start](#) or [My Pathomation](#), that stores slides and provides API which PMA.Studio is utilizing. [More information about what a tile server is can be found on your YouTube channel.](#)

User account

From:

<https://docs.pathomation.com/pimsls/3.0.1/> - **PIMS LS 3.0.1**

Permanent link:

<https://docs.pathomation.com/pimsls/3.0.1/doku.php?id=terminology>

Last update: **2022/11/21 12:23**

