


## Conference settings (administration)

This page allows you to setup PMA.Live server URL that creates live broadcast session for a slide to be viewed by audience using PMA.Studio. You can also validate the PMA.Live URL before saving.



### Configure Conferencing

PMA.Live url

From:

<https://docs.pathomation.com/pimsls/3.1.0/> - **PIMS LS 3.1.0**

Permanent link:

[https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=admin\\_conference](https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=admin_conference)

Last update: **2023/04/06 22:58**

