

# Global settings (administration)


This admin page allows us to set global settings for PMA.Studio application. Page is divided into two set of fields:

**Global settings** has 5 global fields: ·

- 1. Language
- 2. Application IconNo file chosen
- 3. Email For Help
- 4. Conference Server Url
- 5. Snapshot Filename Template

**Annotation settings** has 5 global fields: ·

- 1. Show Annotation Dialog
- 2. Auto Save Annotations
- 3. Restart Annotations Drawing
- 4. Show Slide Info Button
- 5. Show Full Slide Path In Viewport



From:  
<https://docs.pathomation.com/pims/3.1.0/> - **PIMS LS 3.1.0**

Permanent link:  
[https://docs.pathomation.com/pims/3.1.0/doku.php?id=admin\\_global&rev=1648660454](https://docs.pathomation.com/pims/3.1.0/doku.php?id=admin_global&rev=1648660454)

Last update: **2022/03/30 20:14**

