

# Global settings (administration)


This admin page allows us to set global settings for PMA.Studio application. Page is divided into two set of fields:

**Global settings** has 5 global fields: ·

1. Language
2. Application IconNo file chosen
3. Email For Help
4. Conference Server Url
5. Snapshot Filename Template

**Annotation settings** has 5 global fields: ·

1. Show Annotation Dialog
2. Auto Save Annotations
3. Restart Annotations Drawing
4. Show Slide Info Button
5. Show Full Slide Path In Viewport



From:  
<https://docs.pathomation.com/pims/ls/3.1.0/> - **PIMS LS 3.1.0**

Permanent link:  
[https://docs.pathomation.com/pims/ls/3.1.0/doku.php?id=admin\\_global&rev=1648660454](https://docs.pathomation.com/pims/ls/3.1.0/doku.php?id=admin_global&rev=1648660454)

Last update: **2022/03/30 20:14**

