

Global settings (administration)


This admin page allows us to set global settings for PMA.Studio application. Page is divided into two set of fields:

Global settings has 6 fields: ·

1. Language
2. Application IconNo file chosen
3. Email For Help
4. Conference Server Url
5. Snapshot Filename Template

Annotation settings has 5 fields: ·

1. Show Annotation Dialog
2. Auto Save Annotations
3. Restart Annotations Drawing
4. Show Slide Info Button
5. Show Full Slide Path In Viewport



From:
<https://docs.pathomation.com/pimsls/3.1.0/> - PIMS LS 3.1.0

Permanent link:
https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=admin_global&rev=1649165315

Last update: **2022/04/05 16:28**

